HELIOS 700

Wireless Unleashed, Boundless Possibilities

Wireless Intraoral Scanner





Break Free from Cables, Unlock Infinite Possibilities

Say goodbye to cumbersome wired connections, fulfill your scanning needs from all angles, and experience greater ease of use.

High-speed WiFi Connection, Stable Data Transmission

Flagship-level WIFI technology, providing a stable connection with no dropouts.

Zero-latency lossless transmission, faithfully reproducing data.



100_{min}

a single battery can power continuous scanning for (approximately 60+ typical cases)

3100_{mAh}

high-capacity independent dual-battery



Way .

Powerful Battery, Endless Endurance

Superior energy management system, enter sleep mode when docked, instantly awaken when picked up.
Wireless handle charging, "Charge and Go".

Command the Wireless

Unleash the Freedom

High-Precision Reproduction, Seamless Scanning

- · High-performance image sensor combined with Helios exclusive algorithms deliver precision far exceeding clinical requirements.
 - (16µm full-arch accuracy, 8-micron single-tooth accuracy)
- · Upgraded projection chip, 20fps high frame rate for fluid scanning.
- · 24bit digital color innovation, reproducing true intraoral colors.
- · 20mm large depth of field, enhancing imaging capabilities for challenging clinical areas.
- · Automatic adaptive color and precision calibration, eliminating the need for cumbersome calibration tools and processes.

Streamlined Design, Effortless Grip



Exceptional Experience

Light and Extraordinary

Eighteeth Digital Cloud Center

One Cloud, Infinite Insights

Unified Device Ecosystem

Sync data from Helios series, FinScan, and NanoPix series into one intelligent cloud - command your digital workflow from a single cockpit.

Accelerated Amazon Cloud Global Nodes:

Ensuring fast and stable data transmission, supporting remote collaboration and diagnosis worldwide.



Cloud-Local Dual-Mode Orthodontic Simulation:

Multi-scenario orthodontic simulation services that combine the advantages of both cloud and local solutions to meet varied clinical needs.



Helios' orthodontic simulation function can generate online treatment plans, providing clear before-and-after comparisons, helping patients build confidence in the treatment plan.



Model display





Diagnosis details Oral health knowledge

Al personalized oral health reports, helping doctors and patients efficiently communicate diagnoses and treatment plans, bridging the information gap.

The upcoming Eighteeth Digital Cloud Center is coming soon, joining forces with Helios 700 to bring a brand new digital ecosystem experience to dental care.

Diverse Applications, Comprehensively Meeting Clinical Needs



Auto/Manual Bite



Enhanced Implant Workflow



Shining Mode



Occlusion Analysis



Intraoral Camera



Undercut Detection



Texture Mapping



Tooth Measurement



Margin Line Drawing



Software Update



Auto Al



Swap Upper and Lower Jaw



Air Mouse



Model with Base



High-Definition Images

Dual-sided Detachable Scan Tips

- · Safe for 100 high-temperature sterilization cycles.
- · Smooth, edgeless design for a comfortable scanning experience.
- No need for powdering or heating.



Multifunction Buttons

· Dual independent buttons, enabling one-click operation, carefree control.



Silent Fan Anti-fogging

· Silent scanning, optimizing the treatment experience.



Motion Control Free Rotation

· Rotate and inspect scanned data via handle for comprehensive data check.



Recommended System Requirements

Processor	Intel®CoreTM i7-13700H Generation, base frequency 2.4 GHz (or better)
Memory	32GB DDR4, frequency 2666 MHz (or better)
Hard Drive	1T (or more) SSD
Graphics Card	NVIDIA® GeForce® Only (RTX 4060 6G or higher)
Display	15.6"FHD (1920x1080) (or more)
Port	3.0 high-speed port
Operation system	Windows 10/11



Accuracy

Scan Field
Speed
Connection Method
Connectivity Range
Single Battery Capacity
Single Battery Runtime

Full Arch: 16µm Single Tooth: 8µm 16mm×14mm(L) 20fps

Zuips

Dedicated WiFi

6m

3100mAh 100 minutes Number of Batteries
Defogging Method
Calibration

Calibration

Data Output

2 pcs Fan

Precision Calibration: Calibration-free

Calibration-free
Color Calibration:
Calibration-free
STL/PLY/OBJ